

ActarSim Manual

H. Alvarez Pol and D.Y. Pang

May 27, 2009

Abstract

A manual for ActarSim.

Contents

1	Introduction	2
2	Installation	2
2.1	Requirements	2
2.2	Installation	2
2.3	Run the program	3
3	A typical session with ActarSim	3
3.1	Running ActarSim	4
3.2	The digitization macro	6
3.3	Digitization in detail	7
3.3.1	Pads geometry	7
3.3.2	Drift velocity and diffusion parameters	7
3.3.3	Digitization stride-by-stride	8
3.3.4	The wire induction mode	8
3.3.5	A complete example of running the digitizationMacro	9
3.4	An example of analyzing macro	9
3.5	An example of TNtuple reader	11
4	List of messengers in the Geant4 part of ActarSim	11
4.1	Commands controlling the detector description:	12
4.1.1	gas box detector commands	12
4.1.2	silicon detector commands	13
4.1.3	scintillator detector commands	14
4.2	physics processes commands	15

4.3	event generation commands	16
4.3.1	Track a particle or set of particles defined from the Particles list	16
4.3.2	Track a predefined reaction from a file	18
4.3.3	Track a reaction calculated from CINE (program from W. Mittag)	18
4.3.4	Track a reaction calculated from KINE (program from M.S. Golovkov)	19
4.4	event visualization	21
4.5	histogramming control	21
5	Appendix	22
5.1	An example of macro for ActarSim	22

1 Introduction

ActarSim is a code developed for the simulation of the ACTAR experimental setup.

2 Installation

2.1 Requirements

The following software components are required:

GEANT4 - due to the physics definitions, only 4.8.0 and later geant4 versions could be used with this ActarSim simulation. It should be compiled following the instructions provided by the geant4 collaboration (<http://geant4.web.cern.ch/geant4/>). Check the installation with the included geant4 examples.

ROOT - tested on versions 5.06 and some of the latests and the old 4.03. It is mandatory to install the libraries and include directories. Follow the typical installation as shown in <http://root.cern.ch>.

G4UIROOT - not required, but very useful user interface, mainly for beginners because of their fast and painless access to the commands. Download the code from <http://i.home.cern.ch/i/iglez/www/alice/G4UIRoot> and follow the instructions for installation. Check the installation with the modified geant4 examples.

ActarSim - the present code.

2.2 Installation

It is recommended to take the code from the subversion repository using the *subversion* tools. Subversion is available for most of the operative systems, allowing the access to the repository where the code is stored. The required login and password combination could be obtained from the code developers (mail to repository developer, hector.alvarez@usc.es).

In case you have no access to the repository, it is possible to download an stable version of the code Actar-Sim_DATE.tar.gz file on your GEANT working directory. Unpack the file with the commands:

```
tar -xzf ActarSim_DATE.tar.gz
```

A directory ActarSim will be created with the relevant files in.

In any case, after getting the code, one should compile it. To do this, enter in the ActarSim directory and compile it by typing:

```
cd ActarSim
gmake
```

2.3 Run the program

Provided the code was successfully compiled, the code runs using the code name alone (for interactive session) or with a macro name (for a batch session):

```
ActarSim                // interactive session
ActarSim batch1.mac      // batch session
```

Description of the commands which can be used in the macro is given in section 4. It would be useful to specify the terms *number of events* (or *number of Geant4 events*) and *number of reactions* here for the following text. ActarSim calculates information of one *reaction* by issuing two *geant4 events*: one for the beam and a following one for the reaction products (*primaries* in the glossary of Geant4). Since the event numbers in Geant4 start from zero, the **even events correspond to beams and odd events correspond to reaction products**. Having this in mind is important to understand the analyzing macros described in the following sections. Because of this, if a user wants to simulate, say, 100 reactions, he/she has to pass a number 200 to the command `/run/beamOn`:

```
/run/beamOn 200
```

Another thing is that by some reason not known right now, silicon and scintillator information defined in classes *ActarSimSilHit* and *ActarSimSciHit*, respectively, for the i th reaction are written, instead in the $(2i + 1)$ th event, in the $(2i)$ th event (i starts from zero), so if the user really care about the information of the last reaction, for example, the 100th reaction in the previous example, he/she has to pass the number 201 to the command `/run/beamOn`:

```
/run/beamOn 201
```

Of course it does not matter much if we loss information of one reaction if we simulate 10000 reactions.

3 A typical session with ActarSim

A typical ActarSim session contains the following procedures:

1. Run ActarSim, in this step, the
2. Run the digitizationMacro
3. Run the analyzing macro
4. Run the Ntuple reader to analysis the results

3.1 Running ActarSim

An example of macro file needed to run ActarSim is given in section 5.1. Detailed explanations of the commands in this macro file are given in section 4. Most of the commands and their parameters can be kept the same as in this example. The information that has to be provided by the user are the following:

1. The geometry of chamber. Note that the numbers are *half lengths* of the real geometry, for example, for a cubic chamber of $300 \times 300 \times 300 \text{ mm}^3$, the corresponding commands are:

```
/ActarSim/det/gasGeoIncludedFlag on
# if box
/ActarSim/det/gas/setDetectorGeometry box
/ActarSim/det/gas/setXLengthGasBox 150. mm
/ActarSim/det/gas/setYLengthGasBox 150. mm
/ActarSim/det/gas/setZLengthGasBox 150. mm
```

2. To include silicon and scintillator (Csl) detectors or not:

```
/ActarSim/det/silGeoIncludedFlag on or off
/ActarSim/det/sciGeoIncludedFlag on or off
```

At present the geometry of ancillary detectors are hard-coded, i.e., the silicon detectors are $300 \mu\text{m}$ thick and of $100 \times 100 \text{ mm}^2$ square shapes, the square Csl detectors are $25 \times 25 \text{ mm}^2$ and 30 mm thick. If we want to include the silicon and scintillator detectors, we also have to specify the size of “boxes” (usually the same length parameter as the gas chamber) to put the ancillary detectors:

```
/ActarSim/det/sil/sideCoverage 56
/ActarSim/det/sil/xBoxHalfLength 150. mm
/ActarSim/det/sil/yBoxHalfLength 150. mm
/ActarSim/det/sil/zBoxHalfLength 150. mm

/ActarSim/det/sci/sideCoverage 56
/ActarSim/det/sci/xBoxHalfLength 150. mm
/ActarSim/det/sci/yBoxHalfLength 150. mm
/ActarSim/det/sci/zBoxHalfLength 150. mm
```

please refer to section 4 for the definition of *sideCoverage*.

3. The active gas and its pressure, for example, deuterium gas at STP condition:

```
/ActarSim/det/gas/setGasMat D2_STP
```

The following gas/pressure are pre-defined in ActarSim:

- isobutane: isoC4H10STP, isoC4H10_150, isoC4H10_220, isoC4H10_300, isoC4H10_500, isoC4H10_710, isoC4H10_1300, and isoC4H10_1880

- deuterium: D2_40, D2_60, D2_80, D2_100, to D2_400 with steps of 20 mbar, D2_STP, D2_1695, D2_1800, and D2_1950
- helium: He_1900, and He_2010

where the units are in mbar.

4. What information are to be stored in the output file:

```
#Control of the output on the ROOT file
#if all the tracks are stored (note: huge space consumption)
/ActarSim/analControl/storeTracks off
/ActarSim/analControl/storeTrackHistos off
/ActarSim/analControl/storeEvents on
/ActarSim/analControl/storeHistograms off
/ActarSim/analControl/storeSimpleTracks on
```

Usually we do not need to restore all tracks, which consumes huge hard disk space. For the ROOT macros described in the following text, only *storeEvents* and *storeSimpleTracks* are needed to be switched on. With these options, the output file size is around 60 MB for a run of 5000 Geant4 events.

5. Whether we need to treat the beam energy loss (beam interaction) in the gas and its emittance (realistic Beam):

```
/ActarSim/gun/beamInteraction on
/ActarSim/gun/realisticBeam on
/ActarSim/gun/beamRadiusAtEntrance 2.5 mm
/ActarSim/gun/emittance 200.0
```

The unit of emittance here is mm mrad.

6. Which reaction kinematics calculator to use. We have two options: CINE and KINE. Suppose we use KINE, the following commands should be issued:

```
/ActarSim/gun/reactionFromKine on
/ActarSim/gun/Kine/incidentIon 28 78 28 0.0 77.96318
/ActarSim/gun/Kine/targetIon 1 2 1 0.0 2.0141
/ActarSim/gun/Kine/scatteredIon 28 79 28 5.0 78.97107
/ActarSim/gun/Kine/recoilIon 1 1 1 0.0 1.007825
/ActarSim/gun/Kine/labEnergy 624. MeV
/ActarSim/gun/Kine/randomThetaCM on
/ActarSim/gun/Kine/randomThetaRange 0.0 180.0
/ActarSim/gun/Kine/randomPhiAngle on
/ActarSim/gun/Kine/userThetaCM 41.0 deg
/ActarSim/gun/Kine/userPhiAngle 50.0 deg
```

These commands define the entrance- and exit-channel particles to be studied. The parameters for these particles are atomic number, mass number, charge number, excitation energy (in MeV) and atomic mass

(in u), respectively. The above information is for the $^{78}\text{Ni}(d,p)^{79}\text{Ni}$ at incident energy of 8 AMeV with ^{79}Ni at 5 MeV of excitation. The θ_{cm} and ϕ angles are randomized (so that the values of *userThetaCM* and *userPhiAngle* above actually do not matter).

There are also commands that control the visualization of the ActarSim runs, however, they are not relevant to the storage of the calculated informations. Note that enabling the visualization will greatly slow down the speed of calculation.

3.2 The digitization macro

The output file of ActarSim only contain information about strides of particles in the gas (data members of class *ActarSimSimpleTrack*), and about the silicon and CsI detectors (data members of class *ActarSimSilHit* and *ActarSimSciHit*). Pad signals (charges induced and timing, data member of class *ActarPadSignal*) are calculated by using a digitization macro *digitizationMacro.C*. This macro dose the folloing things:

1. reads the output file of ActarSim event-by-event,
2. for each event, reads the information of each stride,
3. for each stride, calculate the corresponding pads that have induced charge on them, for each pad which has induced charge, the corresponding time for the primary electrons of this stride to drift from the ionization position to the pad plane is also calculated.

To perform above calculations, the following information have to be passed to the digitizationMacro:

1. the geometry of the detector:

```
thePadsGeometry.SetGeometryValues(Int_t geometryType,
                                   Int_t padType,
                                   Int_t padLayout,
                                   Double_t xLength,
                                   Double_t yLength,
                                   Double_t zLength,
                                   Double_t radius,
                                   Double_t padSize)
```

Where all distances in this command should be included in mm. The macro calculates the number of pads in each column and row.

2. The drift velocity and diffusion parameters of electron clouds:

```
theDriftManager.SetDriftVelocity(Double_t velocity)      in mm/ns
theDriftManager.SetDiffusionParameters(Double_t long,
                                         Double_t trans) in mm^2/ns
```

3. other parameters, such as the Lorentz angle

```
theDriftManager.SetLorentzAngle(Double_t lor)           in radians
```

The event-loop (and inside a stride-loop) is performed for all strides within all events in the input root file. The projection of each stride on the pad plane is performed using the drift parameters. Then, a calculation of the (relative units) induced charge in each pad close to the projection points is performed:

```
Int_t CalculatePositionAfterDrift(projectionOnPadPlane* pro);
ActarPadSignal* CalculatePadsWithCharge(projectionOnPadPlane* pro,
                                         TClonesArray* clo);
```

Those pads containing induced charge are stored for further analysis in a *TClonesArray* stored in a *TTree* in an output *TFile*. This *TFile* could be used in next analysis steps (the *analysisExample.C*, for example).

3.3 Digitization in detail

3.3.1 Pads geometry

The digitizationMacro calculates the number of pads using an algorithm which divides the available pad plane on the maximum number of pads with the approximate size given by the user. A small change of the pad size could be possible. The configuration and numbering of the pads are the following:

a) square pads: pads are ordered in rows and columns from 1 up to the parameters *numberOfRows* and *numberOfColumns*. An example with 5 columns and 4 rows is shown in Fig.1-(a). In the *thePadsGeometry.SetGeometryValues* function, the *padType* for square pad is 0.

b) hexagonal pads: pads are ordered in rows and columns from 1 up to the parameters *numberOfRows* and *numberOfColumns*. There are two kinds of pad layouts for hexagonal pads: the MAYA-like (Fig.1-(b), *padType*==1 and *padLayout*==0) and the MAYA-tilted layout (Fig.1-(c), *padType*==1 and *padLayout*==1).




Figure 1: Pad types and layout patterns in the simulation code.

In all three cases a single number can identify the pad. The pad number is 1 for row=1, column=1, 2 for row=1, column=2, etc. A set of functions calculates the pad number from the row and column and viceversa.

3.3.2 Drift velocity and diffusion parameters

The electron drift velocity and diffusion parameters can be passed to the digitizationMacro by using the following commands:

```
theDriftManager.SetDriftVelocity(Double_t velocity)      in mm/ns
theDriftManager.SetDiffusionParameters(Double_t long,
                                         Double_t trans) in mm^2/ns
```

For deuterium and isobutane gases, there are analytical formula incorporated in the code to calculate these quantities. To do use these formula, instead of the two commands above, one has to use the following commands instead:

```
theDriftManager.SetDriftParameters(Voltage,height,pressure,gasName)
```

where, for MAYA-like ACTAR detector, *Voltage* (in volts) is the voltage between the upper cathode and the Frish grid, *height* (in mm) is the distance between the upper cathode and the Frish grid, *pressure* (in mbar) is the pressure of the active gas, and *gasName* is the name of the active gas, it has only two optional values: “deuterium” and “isobutane”. At the end of function *theDriftManager.SetDriftParameters()* the drift velocity and diffusion parameters will be automatically calculated.

The analytical formula for deuterium and isobutane gases were obtained by fitting the corresponding curves in ref.[1]. More detailed description of these formula can be found in the *ActarSim-report* by D.Y. Pang.

3.3.3 Digitization stride-by-stride

In digitizationMacro, a loop on all events is made for the Tree in the input TFile. Inside each event, a loop on all *ActarSimSimpleTracks* (strides) is performed. For each stride, an object of the class *projectionOnPadPlane* is assigned. This class contains a pointer to the stride, vectors containing the projection on the pad plane for the initial and final points of the stride and the corresponding drift times for each point. Also the values of the longitudinal and transversal sigma of the diffusion distributions is contained.

All this projections and quantities are obtained from the original stride and the geometrical and drift parameters. Actually is the *driftManager* class which calculates them from the original stride using the function:

```
theDriftManager.CalculatePositionAfterDrift(projection)
```

Next, the charge induced by the electron drift on the pads is calculated. For this, a candidate pad list is prepared including those between the initial and final pads where the projections of the stride lie. A few neighbouring pads are included in the list (by introducing the *securityFactor*) to account for the possible pads diagonal to the initial or final points. Presently, the charge is obtained from the integral of a two dimensional function, which describes the distribution of the induced charges on the pad plane, on the area of the pad (for hexagonal pads, we are using an approximating square limits in the integration – since it is much more complicated to describe the limits of a hexagonal pad). The form of this function depends on the induction mode (wire, MicroMegs, GEM). To perform this calculation, the function:

```
theDriftManager.CalculatePadsWithCharge(projection,padSignalCA)
```

should be called on the event loop.

3.3.4 The wire induction mode

In case of wire amplification, the induced charge distributions, parallel and normal to the wire direction, as a function of the distance between the center of the amplification x , are $\rho_1(x)$ and $\rho_2(x)$, respectively. According to the empirical formula of E. Mathieson and J.S. Gordon [2, 3], $\rho_1(x)$ and $\rho_2(x)$ are:

$$\begin{aligned}\rho(\lambda) &= q_a \times K_1 \frac{1 - \tanh^2(K_2 \lambda)}{1 + K_3 \tanh^2(K_2 \lambda)}, \\ K_1 &= \frac{K_2 \sqrt{K_3}}{4 \tan^{-1} \sqrt{K_3}}, \quad \text{and} \\ K_2 &= \frac{\pi}{2} \left(1 - \frac{\sqrt{K_3}}{2} \right),\end{aligned}\tag{1}$$

where q_a is the net anode charge, $\lambda = x/h$ with h being the distance between the anode wire plane and the cathode pad plane. To calculate the distributions of induced charge on pads using Eq.(1), the only thing needed is the value of Mathieson factor K_3 .

The Mathieson factor K_3 is a function of r_a/s and h/s , where r_a and s are radius of the amplification wire and the distance between two neighbouring wires, respectively. Details of calculations of K_3 can be found in the *ActarSim-report* by D.Y. Pang. If one want to use wire induction mode, in running of the digitizationMacro, one need to issue the following two commands:

```
theAmplificationManager.SetIsWireOn()
theAmplificationManager.SetWireAmplificationParameters(0.02,2.,3.)
```

where the first command switch the wire induction mode on, and the second command passes the wire parameters, r_a , s and h to the program (units in mm). In the above example, $r_a = 20 \text{ } \mu\text{m}$, $s = 2 \text{ mm}$ and $h = 3 \text{ mm}$, respectively.

3.3.5 A complete example of running the digitizationMacro

As explained above, below is one complete example to run the digitizationMacro:

```
$ root -l
root [0] gSystem->Load("actarsim.sl");
root [1] .L digitizationMacro.C++
root [2] thePadsGeometry.SetGeometryValues(0,0,0,150.,150.,150.,100.,2)
root [3] theDriftManager.SetDriftParameters(10000.,300.,1013.25,"deuterium")
root [4] theAmplificationManager.SetIsWireOn()
root [5] theAmplificationManager.SetWireAmplificationParameters(0.02,2.,3.)
root [6] digitEvents("simData.root","digiData.root",0)
```

By doing this, we use the digitizationMacro to read the ActarSim output file *simData.root* and write the resulting digitization file into *digiData.root*. The active target gas is deuterium at 1 atm pressure. We are using wire induction mode with wire radius of 0.02 mm, intervals between wires are 2 mm and the wire plane is 3 mm above the pad plane.

3.4 An example of analyzing macro

After the digitization, we have the *simulation data* which is equivalent to the experimental data. The simulation data is stored in two file: ancillary detector signals (silicon and CsI) are stored in the output file of ActarSim (*simData.root*) and pad signals (induced charge and time) are stored in the output file of the digitizationMacro (*digiData.root*). Both files are written in a event-by-event way and events related to the same reaction have the same event ID in both files.

In principle, all algorithms developed in the analysis of the MAYA experiment should be able to be applied to the simulation data in the *simData.root* and *digiData.root* files. We provide an example macro *analysisExample.C* for this purpose.

For ^{78}Ni induced (d, p) , (d, d') and (d, t) reactions, the analysisExample macro reconstruct the trajectory of the light particle, i.e., proton, deuteron, and triton. From each reaction that the trajectory of the light particle is reconstructable, the analysisExample macro reads the energy signal in silicon and CsI detectors and calculates the range of that particle inside the gas. Once the trajectory is known, the θ_{Lab} and ϕ angle of the out-going

particle in the laboratory system are able to be calculated and compared with their corresponding “real” values and then we are able to analyse the angular resolution of the detector for a given gas/pressure, pad information (pad size, shape and pattern of layout), and induction mode. Similarly, the position resolution can also be obtained by comparing the reconstructed reaction vertex Z-value with its “real” value given by Geant4.

To run this macro, the following ROOT command should be used:

```
$ root -l
root [0] gSystem->Load("actarsim.sl");
root [1] .L analysisExample.C++
root [2] thePadsGeometry.SetGeometryValues(0,0,0,150.,150.,150.,100.,2)
root [3] theDriftManager.SetDriftParameters(10000.,300.,1013.25,"deuterium")
root [4] reader("simData.root","digiData.root","Ntuples.root",1000.,3,2,0,0,0)
```

Note that the arguments for functions *SetGeometryValues()* and *SetDriftParameters()* should be the same as the ones used in the corresponding digitizationMacro. The arguments for the function *reader()* are

1. *simData.root*: the output file of ActarSim
2. *digiData.root*: the output file of the digitizationMacro
3. *Ntuples.root*: the output of this analysisExample macro, in which the following information are stored reaction-by-reaction in a TNtuple:
 - G4VertexZ: the Z-value of the reaction vertex given by Geant4
 - G4ThetaLab: the θ_{Lab} of the out-going light particle given by Geant4
 - G4PhiLab: the ϕ angle of the out-going light particle given by Geant4
 - G4VertexE: the beam energy at the reaction vertex given by Geant4 (taking into account the beam energy loss in the gas)
 - calVertexZ: the Z-value of the reaction vertex from the reconstructed particle trajectory
 - calThetaLab: the θ_{Lab} of the out-going light particle from the reconstructed particle trajectory
 - calPhiLab: the ϕ angle of the out-going light particle from the reconstructed particle trajectory
 - eSil: the energy loss of the light particle inside the silicon detector (in MeV)
 - eSci: the energy loss of the light particle inside the CsI detector (in MeV)
 - protonRangeInGas: the proton range inside the gas (in mm), note that the range is calculated whenever the trajectory of the particle is reconstructable, it does not mean that the particle is stopped inside the gas.
4. the dynamic range: ratio between the maximum and minimum charge signals that can be measured on the pads. In this analysisExample macro, we determine the minimum charge on the pads by dividing the maximum charge with the dynamic range. This is equivalent to the thresholds. In the above example, dynamic range is 1000.
5. the beam projection width (bw): half of the number of rows of pads that are taken by the beam, the trajectory of light particles is not able to be reconstructed if its projection on pad plane falls between $[\text{numberOfRows}/2 - \text{bw}, \text{numberOfRows}/2 + \text{bw}]$. In the above example, bw is 3.
6. the margin width (mw): number of rows and columns at the edges of the gas chamber. We do not take the signals of the “margin pads” as effective signals. In the above example, mw is 2.

7. the last three zeroes in the arguments of the function `reader()` are verbose level, minimum and maximum reaction numbers (`minReactionNumber` and `maxReactionNumber`), respectively. If the verbose level is larger or equal to 1, the program will give some diagnosis information, which is useful if we want to know the details of the running of this macro. If `minReactionNumber` and `maxReactionNumber` are zeroes, all reactions will be treated, otherwise, only reactions with reaction number between the `minReactionNumber` and `maxReactionNumber` will be treated. In the above example, we include all reactions in the analysis.

In this `analysisExample` macro, seven histograms are defined, namely:

- `hLightEvsTheta`: two dimensional, $E - \theta_{\text{Lab}}$ kinematics plot for the out-going light particle, where E are energy signals in ancillary detectors.
- `hDThetaVsTheta`: two dimensional, $\Delta\theta_{\text{Lab}} - \theta_{\text{Lab}}$, this histogram shows how the difference between “real” value of θ_{Lab} and its reconstructed value ($\Delta\theta_{\text{Lab}}$) depends on the value θ_{Lab} itself.
- `hDThetaVsPhi`: two dimensional, $\Delta\theta_{\text{Lab}} - \phi$, this histogram shows how $\Delta\theta_{\text{Lab}}$ depends on the reconstructed ϕ value.
- `hRangeVSTheta`: two dimensional, correlation between the range of the light particle and its θ_{Lab} .
- `hDTheta`: one dimensional, angular resolution in θ_{Lab}
- `hDPhi`: one dimensional, angular resolution in ϕ
- `hDVertexZ`: one dimensional, position resolution in reaction vertex Z -value.

The users of ActarSim can define their own histograms. It is more useful to study the correlations between quantities we are interested in. This can be done by using `TNtuple`s.

3.5 An example of `TNtuple` reader

The `analysisExample` macro generate an output ROOT file which contains one `TNtuple`. The contents of this `ntuple` has been introduced in the previous section. We provide the macro `ntupleReader.C` as an example of manipulate these `ntuple` data members. In this macro we can include as many number of runs of simulation and analysis them together. For example, the (d, p) , (d, d') and (d, t) reactions are simulated in different runs, we can use this macro to put analysis these reactions together.

Some example of the analyzing result using the `analysisExample` and `ntupleReader` macros can be find in the ⁷⁸*Ni-dp-simulation-report* file by D.Y. Pang.

4 List of messengers in the Geant4 part of ActarSim

Next, a complete description of the possible messenger commands available for selecting different options and for the modification of the behavior of the program. The different messenger commands are ordered and separated by their utility, in different directories. The directories are:

- **det**: containing the detector commands, able to modify the setup conditions. It contains two subdirectories: `det/sil` controlling the Silicon Detectors setup and `det/sci` controlling the Scintillator Detectors setup.
- **phys**: containing all related with the interaction physics;

- **gun:** commands related with the event generation;
- **event:** event visualization;
- **analControl:** histogramming control.

For each case, the command complete name initiates the description, following by the command explanation and some additional features, if required.

4.1 Commands controlling the detector description:

```
/ActarSim/det/
    Command description.  Just the directory name...

/ActarSim/det/gasGeoIncludedFlag
    Includes the geometry of the gas volume in the simulation (default "off").

/ActarSim/det/silGeoIncludedFlag
    Includes the geometry of the silicons in the simulation (default "off").

/ActarSim/det/sciGeoIncludedFlag
    Includes the geometry of the scintillator in the simulation (default "off").

/ActarSim/det/setMediumMat
    Select Material of the Medium.

/ActarSim/det/setEleField
    Define electric field.
    Usage: /ActarSim/det/setEleField Ex Ey Ez (in MV/mm)

/ActarSim/det/setMagField
    Define magnetic field.
    Usage: /ActarSim/det/setMagField Bx By Bz unit

/ActarSim/det/update
    Update geometry.
    This command MUST be applied before \"beamOn\"
    if you changed geometrical value(s) or gas density.

/ActarSim/det/print
    Prints geometry.
```

4.1.1 gas box detector commands

```
/ActarSim/det/setGasMat
    Select Material of the Gas.
    DefaultValue("isoC4H10")
```

```

/ActarSim/det/gas/setBeamShieldMat
    Select Material of the beam shield.
    DefaultValue("isoC4H10")

/ActarSim/det/gas/setDetectorGeometry
    Select the geometry of the detector.
    DefaultValue("box")

/ActarSim/det/gas/setBeamShield
    Sets a beam shield and selects the geometry.
    DefaultValue("tube")

/ActarSim/det/gas/setBeamShieldMat
    Select Material of the beam shield.
    DefaultValue("isoC4H10")

/ActarSim/det/gas/setXLengthGasBox
    Select the half-length X dimension of the Gas Box

/ActarSim/det/gas/setYLengthGasBox
    Select the half-length Y dimension of the Gas Box

/ActarSim/det/gas/setZLengthGasBox
    Select the half-length Z dimension of the Gas Box

/ActarSim/det/gas/setRadiusGasTub
    Select the external radius of the Gas Tube.

/ActarSim/det/gas/setLengthGasTub
    Select the half-length of the Gas Tube.

/ActarSim/det/gas/setInnerRadiusBeamShieldTub
    Select the external radius of the Gas Tube.

/ActarSim/det/gas/setRadiusBeamShieldTub
    Select the internal radius of the Gas Tube.

/ActarSim/det/gas/setLengthBeamShieldTub
    Select the half-length of the Gas Tube.

/ActarSim/det/gas/print
    Prints geometry.

```

4.1.2 silicon detector commands

```

/ActarSim/det/sil/
    Silicon detector control

```

/ActarSim/det/sil/print
Prints geometry.

/ActarSim/det/sil/sideCoverage
Selects the silicon coverage (default 1).
6 bits to indicate which sci wall is present (1) or absent (0).
The order is:
bit1 (lsb) beam output wall
bit2 lower (gravity based) wall
bit3 upper (gravity based) wall
bit4 left (from beam point of view) wall
bit5 right (from beam point of view) wall
bit6 (msb) beam entrance wall
Convert the final binary to a decimal number and set the coverage!

/ActarSim/det/sil/xBoxHalfLength
Sets the x half length of the silicon detectors box
DefaultValue(0.5)

/ActarSim/det/sil/yBoxHalfLength
Sets the y half length of the silicon detectors box
DefaultValue(0.5)

/ActarSim/det/sil/zBoxHalfLength
Sets the z half length of the silicon detectors box
DefaultValue(0.5)

4.1.3 scintillator detector commands

/ActarSim/det/sci/
Scintillator detector control

/ActarSim/det/sci/print
Prints geometry.

/ActarSim/det/sci/sideCoverage
Selects the scintillator coverage (default 1).
6 bits to indicate which sci wall is present (1) or absent (0).
The order is:
bit1 (lsb) beam output wall
bit2 lower (gravity based) wall
bit3 upper (gravity based) wall
bit4 left (from beam point of view) wall
bit5 right (from beam point of view) wall
bit6 (msb) beam entrance wall
Convert the final binary to a decimal number and set the coverage!

```
/ActarSim/det/sci/xBoxHalfLength
    Sets the x half length of the scintillator detectors box
    DefaultValue(0.5)
```

```
/ActarSim/det/sci/yBoxHalfLength
    Sets the y half length of the scintillator detectors box
    DefaultValue(0.5)
```

```
/ActarSim/det/sci/zBoxHalfLength
    Sets the z half length of the scintillator detectors box
    DefaultValue(0.5)
```

IMPORTANT NOTE: */ActarSim/det/update* should be called after each geometrical modification. This command **MUST** be applied before */run/beamOn*.

Note: the geometries of sil and sci detectors are hard-coded right now, it should be modified to cope with different setups and possibilities.

4.2 physics processes commands

The full bank of particles (*ConstructBosons()*; *ConstructLeptons()*; *ConstructMesons()*; *ConstructBaryons()*; *G4ShortLivedConstructor* and *G4IonConstructor*) are constructed. The physics processes (transportation, em and decay) are now taken from the *examples/extended/medical/GammaTherapy* (basically the *PhysicsList* implementation for em processes). Other hadronic particles and processes are also taken from this example. Note that the *PhysicsList* can be selected between several possibilities:

```
    em: standard, lowenergy, penelope, (choose one from these three)
    common: decay
    hadronic: elastic, binary, binary_ion, gamma_nuc
    ion low-energy: ion-LowE, ion-LowE-ziegler1977, ion-LowE-ziegler1985,
                  ion-LowE-ziegler2000, ion-standard
```

```
/ActarSim/phys/
    Physics list commands directory
```

```
/ActarSim/phys/setGCut
    Set gamma cut.
```

```
/ActarSim/phys/setECut
    Set electron cut.
```

```
/ActarSim/phys/setPCut
    Set positron cut.
```

```
/ActarSim/phys/setCuts
    Set cut for all.
```

```
/ActarSim/phys/addPhysics
    Add modula physics list.
```

```
/ActarSim/phys/verbose
    Set verbose level for processes
```

Note that a */run/initialization* is now needed to initialize the application. This initialization has been removed from the *ActarSim.cc* main function to permit the definition of different physics (lowenergy, penelope...) from the PhysicsList before the initialization. An example of how to initialize is shown in the macros (test1.mac, for instance).

4.3 event generation commands

The event generation can be controled by user commands under the directory */ActarSim/gun* and subdirectories inside. There are several types of reactions that the program handles:

A) Track a particle or set of particles defined from the Particles list.

This mode is selected setting the commands

```
/ActarSim/gun/reactionFromFile off (default behavior)
/ActarSim/gun/reactionFromCine off
...
```

B) Track a predefined reaction from a file

C) Track a reaction calculated from CINE (program from W. Mittig)

D) Track a reaction calculated from KINE (program from M.S. Golovkov)

For each case, a realistic interaction (beam distribution and interaction over all gas volume, plus the energy lost of the beam in the gas before the reaction). These commands are described below:

```
/ActarSim/gun/realisticBeam
    Simulates beam emittance according to emittance parameters.
    Choice : on, off(default)
```

```
/ActarSim/gun/beamInteraction
    Simulates the beam energy loss in gas.
    Choice : on, off(default)
```

```
/ActarSim/gun/emittance
    Selects the value of the emittance [in mm mrad].
    Default value is 1 mm mrad.
```

```
/ActarSim/gun/beamRadiusAtEntrance
    Selects the beam radius at entrance of ACTAR.
    Used with the emittance to calculate the position and angle
    distributions of the beam when a realisticBeam option is set.
```

4.3.1 Track a particle or set of particles defined from the Particles list

```
/ActarSim/gun/List
```


List available particles. Invoke G4ParticleTable.

/ActarSim/gun/particle

Select the incident particle (proton is default)
(ion can be specified for shooting ions).

/ActarSim/gun/ion

Set properties of ion to be generated.
[usage] /gun/ion Z A Q E
Z:(int) AtomicNumber
A:(int) AtomicMass
Q:(int) Charge of Ion (in unit of e)
E:(double) Excitation energy (in keV)

/ActarSim/gun/energy

Sets the kinetic energy of the primary particle
DefaultValue(20.)

/ActarSim/gun/direction

Set momentum direction.
Direction needs not to be a unit vector.

/ActarSim/gun/position

Set starting position of the particle.

/ActarSim/gun/time

Set initial time of the particle.

/ActarSim/gun/polarization

Set polarization.

/ActarSim/gun/number

Set number of particles to be generated.

/ActarSim/gun/randomVertexZPosition on

Set the vertex Z-value random or not
Choice: on (default), off

/ActarSim/gun/randomVertexZRange min max unit

This command is effective only if randomVertexZPosition=="on"
min: minimum vertex Z-value, default: 0
max: maximum vertex Z-value, default: 300
unit: default mm

/ActarSim/gun/vertexZPosition z-value unit

This command is effective only if randomVertexZPosition=="off"
z-value: user-given vertex Z-value

unit: default mm

4.3.2 Track a predefined reaction from a file

```
/ActarSim/gun/reactionFromFile
    Select a reaction from an input file
    Choice : on, off(default)
```

```
/ActarSim/gun/reactionFile
    Select the reaction definition file
```

4.3.3 Track a reaction calculated from CINE (program from W. Mittig)

```
/ActarSim/gun/reactionFromCine
    Select a reaction using Cine
    Choice : on, off(default)
```

Note that the following commands are under the subdirectory Cine

```
/ActarSim/gun/Cine/randomTheta
    Select a random Theta angle for the scattered particle.
    Choice : on(default), off
```

```
/ActarSim/gun/Cine/randomThetaVal
    Sets the limit in the Theta angle for the scattered particle.
    The value is randomly chosen between the limits (degrees)
    [usage] /ActarSim/gun/Cine/incidentIon min max
    min: the minimum angle (default 0 degrees)
    max: the maximum angle (default 180 degrees)
```

```
/ActarSim/gun/Cine/thetaLabAngle
    Sets theta lab angle for the scattered particle (degrees)
    DefaultValue(0.5)
```

```
/ActarSim/gun/Cine/incidentIon
    Set properties of incident ion to be generated.
    [usage] /ActarSim/gun/Cine/incidentIon Z A Q E
    Z:(int) AtomicNumber
    A:(int) AtomicMass
    Q:(int) Charge of ion (in unit of e)
    E:(double) Excitation energy (in keV)
```

```
/ActarSim/gun/Cine/targetIon
    Set properties of target ion to be generated.
    [usage] /ActarSim/gun/Cine/targetIon Z A Q E
    Z:(int) AtomicNumber
    A:(int) AtomicMass
```

Q:(int) Charge of ion (in unit of e)
E:(double) Excitation energy (in keV)

/ActarSim/gun/Cine/scatteredIon

Set properties of scattered ion to be generated.

[usage] /ActarSim/gun/Cine/scatteredIon Z A Q E
Z:(int) AtomicNumber
A:(int) AtomicMass
Q:(int) Charge of ion (in unit of e)
E:(double) Excitation energy (in keV)

/ActarSim/gun/Cine/recoilIon

Set properties of recoil ion to be generated.

[usage] /ActarSim/gun/Cine/recoilIon Z A Q E
Z:(int) AtomicNumber
A:(int) AtomicMass
Q:(int) Charge of ion (in unit of e)
E:(double) Excitation energy (in keV)

/ActarSim/gun/Cine/reactionQ

Sets the reaction Q-value for the ground state (MeV)

/ActarSim/gun/Cine/labEnergy

Sets the laboratory energy (MeV)

4.3.4 Track a reaction calculated from KINE (program from M.S. Golovkov)

/ActarSim/gun/reactionFromKine

Select a reaction using KINE

Choice : on(default), off

Note that the following commands are under the subdirectory Kine

/ActarSim/gun/Kine/randomTheta

Select a random Theta angle for the scattered particle.

Choice : on(default), off

/ActarSim/gun/Kine/randomThetaVal

Set the range of theta angles if it is to be randomly chosen

[usage] /ActarSim/gun/Kine/randomThetaVal min, max
min: the minimum angle (default 0 degrees)
max: the maximum angle (default 180 degrees)

/ActarSim/gun/Kine/incidentIon

Set the properties of the incident ion

[usage] /ActarSim/gun/Kine/incidentIon Z A Q Ex Mass
Z: atomic number (I)

A: mass number (I)
Q: charge number (I)
Ex: excitation energy (D, in MeV)
Mass: mass (D, in u)

/ActarSim/gun/Kine/targetIon

Set the properties of the target ion

[usage] /ActarSim/gun/Kine/targetIon Z A Q Ex Mass
Z: atomic number (I)
A: mass number (I)
Q: charge number (I)
Ex: excitation energy (D, in MeV)
Mass: mass (D, in u)

/ActarSim/gun/Kine/scatteredIon

Set the properties of the scattered ion

[usage] /ActarSim/gun/Kine/scatteredIon Z A Q Ex Mass
Z: atomic number (I)
A: mass number (I)
Q: charge number (I)
Ex: excitation energy (D, in MeV)
Mass: mass (D, in u)

/ActarSim/gun/Kine/recoilIon

Set the properties of the recoiled ion

[usage] /ActarSim/gun/Kine/recoilIon Z A Q Ex Mass
Z: atomic number (I)
A: mass number (I)
Q: charge number (I)
Ex: excitation energy (D, in MeV)
Mass: mass (D, in u)

/ActarSim/gun/Kine/labEnergy

Set the incident energy

[usage] /ActarSim/gun/Kine/labEnergy E MeV
E: incident energy (D)
MeV: the unit

/ActarSim/gun/Kine/randomThetaCM

Randomize or not the Theta_cm value

Choice: on (default), off

Default: on

/ActarSim/gun/Kine/randomThetaRange min max unit

This command is effective only if randomThetaCM=="on"

min: default value 0.

max: default value 180.

```

    unit: default deg

/ActarSim/gun/Kine/randomPhiAngle
    Randomize or not the Phi angle, useful when testing analyzing macros
    Choice: on (default), off
    default: on

/ActarSim/gun/Kine/userThetaCM theta_cm unit
    This command is effective only if randomThetaCM is "off"
    theta_cm: user given Theta_cm for KINE
    unit: default deg

/ActarSim/gun/Kine/userPhiAngle Phi unit
    This command is effective only if randomPhiAngle is "off"
    Phi: user given phi angle of the recoiled particle,
        e.g., proton in the  $^{78}\text{Ni}(d,p)^{79}\text{Ni}$  reaction in inverse kinematics
    unit: default deg

```

4.4 event visualization

Inside the *ActarSimEventAction* one can chose the visualization of the tracks in the event, in particular the user command:

```

/ActarSim/event/drawTracks
    Draw the tracks in the event (Choice : none, charged, neutral, all(default))
    DefaultValue("all");

/ActarSim/event/printModulo
    Print events (modulo n, that is prints every n event)

```

4.5 histogramming control

The Tree and histograms can be controlled by user commands under the directory */ActarSim/analControl* and subdirectories inside:

```

/ActarSim/analControl/storeTracks
    Store the tracks in the output Tree
    Choice : on, off(default)

/ActarSim/analControl/storeTrackHistos
    Store the tracks in Histograms
    Choice : on, off(default)

/ActarSim/analControl/storeEvents
    Store the events in the output Tree
    DefaultValue("off")

```

```
/ActarSim/analControl/storeSimpleTracks
    Store the simple tracks in the output Tree
    DefaultValue("off")
```

```
/ActarSim/analControl/storeHistograms
    Store the events in the output Tree
    Choice : on, off(default)
```

5 Appendix

5.1 An example of macro for ActarSim

This macro simulates 5000 $^{78}\text{Ni}(\text{d},\text{p})^{79}\text{Ni}$ reactions at 8 A MeV with ^{79}Ni at an excited state of 5 MeV.

```
#####
#*-- AUTHOR : Hector Alvarez-Pol
#*-- Date: 05/2005
#*-- Last Update: 15/05/08
#*-- Copyright: GENP (Univ. Santiago de Compostela)
# -----
# Comments:
#   - 15/05/08 Multidetector geometry
#   - 05/05/06 Updating to new ActarSim (geant4.8) code
#   - 22/11/05 Updated including CINE options
#   - 31/05/05 Macro for ACTAR simulation
#
#####
# Macro file for testing online jobs
#####
# verbosity levels and saveHistory
/control/verbose 0
/control/saveHistory
/run/verbose 0
/event/verbose 0
/tracking/verbose 0
#
# Setting the Physics modules; valid values are here listed:
#   em: standard, lowenergy, penelope, (choose one from this three)
#   common: decay,
#   hadronic: elastic, binary, binary_ion, gamma_nuc,
#   ion low-energy: ion-LowE, ion-LowE-ziegler1977, ion-LowE-ziegler1985,
#   ion-LowE-ziegler2000, ion-standard
#
/ActarSim/phys/addPhysics standard
#/ActarSim/phys/addPhysics decay
#/ActarSim/phys/addPhysics elastic
```

```

#/ActarSim/phys/addPhysics binary
#/ActarSim/phys/addPhysics binary_ion
#/ActarSim/phys/addPhysics gamma_nuc
#/ActarSim/phys/addPhysics lowenergy
#/ActarSim/phys/addPhysics ion-LowE
#/ActarSim/phys/addPhysics ion-LowE-ziegler1977
#/ActarSim/phys/addPhysics ion-LowE-ziegler1985
#/ActarSim/phys/addPhysics ion-LowE-ziegler2000
#/ActarSim/phys/addPhysics ion-standard
#/ActarSim/phys/addPhysics penelope
#
# Cuts for the particles (incomplete list, see README)
#
# /ActarSim/phys/setGCut 0.1
# /ActarSim/phys/setECut 0.1
# /ActarSim/phys/setPCut 0.1
# /ActarSim/phys/setCuts 0.1
/ActarSim/phys/verbose 0
#
# Initialization is moved here from the main allowing PhysicsList
#
/run/initialize
#
# DETECTOR CHARACTERIZATION
#
# GENERAL COMMANDS
#
# Control of the materials
/ActarSim/det/setMediumMat Galactic
#Electric and Magnetic fields
/ActarSim/det/setEleField 0 0 0
/ActarSim/det/setMagField 0 0 0 T
#
# GAS DETECTOR
#
/ActarSim/det/gasGeoIncludedFlag on
# if box
/ActarSim/det/gas/setDetectorGeometry box
/ActarSim/det/gas/setXLengthGasBox 150. mm
/ActarSim/det/gas/setYLengthGasBox 150. mm
/ActarSim/det/gas/setZLengthGasBox 150. mm
#
# if tube
# /ActarSim/det/gas/setDetectorGeometry tube
# /ActarSim/det/gas/setRadiusGasTub 1.0 cm
# /ActarSim/det/gas/setLengthGasTub 1.0 cm
#

```

```

# Beam shield? tube or off
/ActarSim/det/gas/setBeamShield off
# /ActarSim/det/gas/setBeamShield on
# /ActarSim/det/gas/setBeamShieldMat Ion
# /ActarSim/det/gas/setInnerRadiusBeamShieldTub 1.0 mm
# /ActarSim/det/gas/setRadiusBeamShieldTub 1.001 mm
# /ActarSim/det/gas/setLengthBeamShieldTub 1.0 m
#
# gas material: isoC4H10STP, isoC4H10_150, isoC4H10_220,
#               isoC4H10_300, isoC4H10_500, isoC4H10_710,
#               isoC4H10_1300, isoC4H10_1880
#               D2_40, D2_60, D2_80, D2_100, to D2_400 with steps of 20 mbar
#               D2_STP, D2_1695, D2_1800, D2_1950
#               He_1900, He_2010
#
/ActarSim/det/gas/setGasMat D2_STP
#
# SILICON DETECTOR
#
/ActarSim/det/silGeoIncludedFlag on
#Options for Silicon and scintillator coverage:
# 6 bits to indicate which sci wall is present (1) or absent (0)
# order is:
# bit1 (lsb) beam output wall          1
# bit2 lower (gravity based) wall      2
# bit3 upper (gravity based) wall      4
# bit4 left (from beam point of view) wall  8
# bit5 right (from beam point of view) wall 16
# bit6 (msb) beam entrance wall       32
# examples: 63 full coverage; 3 only output and bottom walls ...
/ActarSim/det/sil/sideCoverage 56
/ActarSim/det/sil/xBoxHalfLength 150. mm
/ActarSim/det/sil/yBoxHalfLength 150. mm
/ActarSim/det/sil/zBoxHalfLength 150. mm
/ActarSim/det/sil/print
#
# SCINTILLATOR DETECTOR
#
/ActarSim/det/sciGeoIncludedFlag on
# see above explanation in the equivalent command for the Silicons
/ActarSim/det/sci/sideCoverage 56
/ActarSim/det/sci/xBoxHalfLength 150. mm
/ActarSim/det/sci/yBoxHalfLength 150. mm
/ActarSim/det/sci/zBoxHalfLength 150. mm
/ActarSim/det/sci/print
#
#Control of the output on the ROOT file

```



```

#all the tracks are stored (note: huge space consumption)
#Note: it should come before the update!!!
/ActarSim/analControl/storeTracks off
/ActarSim/analControl/storeTrackHistos off
/ActarSim/analControl/storeEvents on
/ActarSim/analControl/storeHistograms off
/ActarSim/analControl/storeSimpleTracks on
#/ActarSim/analControl/setMinStrideLength 0.1
#
# Update is mandatory after any material,field or detector change
#
/ActarSim/det/update
/ActarSim/det/print
#
# Control of the primary events
#For all cases the possibility to have realistic beam distribution
/ActarSim/gun/beamInteraction on
/ActarSim/gun/realisticBeam on
/ActarSim/gun/beamRadiusAtEntrance 2.5 mm
/ActarSim/gun/emittance 200.0
#
# Realistic Event-Generator on
/ActarSim/gun/reactionFromEvGen off
#
# Reaction from Cine and Event-Generator
/ActarSim/gun/reactionFromCine off
/ActarSim/gun/reactionFromCrossSection off
#
# A) Track a particle or set of particles defined from the Particles list
#
#/ActarSim/gun/List
#/ActarSim/gun/particle proton
#if you want to use an ion, write "ion" in the previous command
#and set the Z, A and charge state in the next...
#/ActarSim/gun/ion 3 11 3
#/ActarSim/gun/energy 11. MeV
#/ActarSim/gun/direction 0 0 1
#/ActarSim/gun/time 0
#/ActarSim/gun/polarization 0
#/ActarSim/gun/number 1
/ActarSim/gun/randomVertexZPosition on
/ActarSim/gun/randomVertexZRange 145. 155. mm
/ActarSim/gun/vertexZPosition 2.0 cm
#
#
# B) Track a predefined reaction from a file:
#

```

```

/ActarSim/gun/reactionFromFile off
#if you select a reaction, you should say the file in the command below
#possibilities are (up to now): He8onp_100MeV_Elastic.dat,
#He8onp_100MeV_tritium.dat, He8onC12_100MeV_Elastic.dat,
#He8onp_50MeV_Elastic.dat, He8onC12_50MeV_Elastic.dat, He8onp_50MeV_tritium.dat
#
#/ActarSim/gun/reactionFile /data/He8onp_50MeV_Elastic.dat
#/ActarSim/gun/randomReaction off      STILL NOT DONE -- ALWAYS RANDOM
#/ActarSim/gun/rowInFileReaction 4     STILL NOT DONE -- ALWAYS RANDOM
#
#C) Track a reaction calculated from CINE (program from W. Mittig)
#
/ActarSim/gun/reactionFromCine off
#/ActarSim/gun/Cine/incidentIon 3 8 3 0.0
#/ActarSim/gun/Cine/targetIon 2 4 2 0.0
#/ActarSim/gun/Cine/scatteredIon 3 8 3 1.0
#/ActarSim/gun/Cine/recoilIon 2 4 2 0.0
#/ActarSim/gun/Cine/reactionQ 1. MeV
#/ActarSim/gun/Cine/labEnergy 70 MeV
#/ActarSim/gun/Cine/randomTheta on
#/ActarSim/gun/Cine/thetaLabAngle 0.1
#
#D) Track a reaction calculated from KINE
#
/ActarSim/gun/reactionFromKine on
/ActarSim/gun/Kine/incidentIon 28 78 28 0.0 77.96318
/ActarSim/gun/Kine/targetIon 1 2 1 0.0 2.0141
/ActarSim/gun/Kine/scatteredIon 28 79 28 5.0 78.97107
/ActarSim/gun/Kine/recoilIon 1 1 1 0.0 1.007825
/ActarSim/gun/Kine/labEnergy 624. MeV
/ActarSim/gun/Kine/randomThetaCM on
/ActarSim/gun/Kine/randomThetaRange 0.0 180.0
/ActarSim/gun/Kine/randomPhiAngle on
/ActarSim/gun/Kine/userThetaCM 41.0 deg
/ActarSim/gun/Kine/userPhiAngle 50.0 deg
#
# VISUALIZATION
#
# Draw the whole geometry tree with details as function of verbosity
# /vis/ASCIITree/verbose 10
# /vis/drawTree
# visualization
# /vis/scene/create
# /vis/open OGLIX
# /vis/viewer/set/autoRefresh 0
#/vis/viewer/flush
# set camera

```

```

# /vis/viewer/reset
# /vis/viewer/set/hiddenEdge 1
# /vis/viewer/set/lightsThetaPhi 120 40
#/vis/viewer/set/viewpointThetaPhi 115. 145.
# /vis/viewer/set/viewpointThetaPhi 90. 90.
# /vis/viewer/zoom 1.0
# /vis/viewer/set/background 1 1 1 1
#/vis/viewer/flush
#
# drawing style
# /vis/viewer/set/style surface
#/vis/viewer/set/style wireframe
#/vis/viewer/flush
#
# drawing the tracks
# /tracking/storeTrajectory 1
# /vis/scene/endOfEventAction accumulate
# /vis/viewer/set/autoRefresh 1
#
# create an empty scene and add the detector geometry to it
#/vis/drawVolume
#/vis/scene/add/axes 0 0 0 0.1 m
#/ActarSim/event/drawTracks all
#/ActarSim/event/printModulo 1
#
# RUN: number of events
#
/run/beamOn 10002

```

Index

ActarSim session, [3](#)
ancillary detector, [9](#)
ancillary detectors, [4](#)
angular resolution, [10](#), [11](#)

beam energy loss, [5](#)
beam projection width, [10](#)

CINE, [5](#), [18](#)

detector commands, [11–14](#)
deuterium, [5](#)
diffusion parameters, [6](#), [7](#)
digitization, [6](#)
drift velocity, [6](#), [7](#)
dynamic range, [10](#)

emittance, [5](#)
event generation, [12](#), [16](#)

gas/pressure, [4](#)
geometry, [4](#)

helium, [5](#)
hexagonal pads, [7](#)
histogram, [11](#)
histogramming control, [12](#), [21](#)

installation, [2](#)
interaction physics, [11](#), [15](#)
isobutane, [4](#)

KINE, [5](#), [19](#)

margin width, [10](#)
Mathieson factor, [9](#)
messenger commands, [11](#)

ntuple, [3](#), [10](#), [11](#)
number of events, [3](#)
number of reactions, [3](#)

pad layouts, [7](#)
pad signals, [9](#)
position resolution, [10](#), [11](#)

range, [9](#)

square pads, [7](#)

stride, [7](#)
subversion repository, [2](#)

the digitizationMacro, [6](#), [9](#)
thresholds, [10](#)
trajectory, [9](#)

visualization, [12](#), [21](#)

wire induction mode, [8](#)

References

- [1] Anna Peisert and Fabio Sauli, Drift and Diffusion of electrons in gases: a compilation, CERN 84-08, 1984
- [2] E. Mathieson and J.S. Gordon, Nucl. Instru. Meth. 227 (1984) 277.
- [3] E. Mathieson, Nucl. Instru. Meth. **A270** (1988) 602.